

# IMPLANT LINE-UP

## **TRAUS SIP20**

### **Premium Implant Engine**



SIP20

**EZ OPEN**

**80 N·cm**

**OPTIC**

- Maintenance friendly
- Intuitive user interface
- Real time RPM display
- Motor auto-calibration
- Efficient gear system

|      |        |                |
|------|--------|----------------|
| 20:1 | RPM    | 15 ~ 2,000 rpm |
|      | Torque | 5 ~ 80 N·cm    |



## **TRAUS SIP10**

### **Advanced Implant Engine**



SET UP


**COMPACT**

**70 N·cm**

**OPTIC**

- 9 programmable memories
- Smart error diagnosis system
- Real time torque & RPM display
- Motor auto-calibration
- Ergonomic design

|      |        |                |
|------|--------|----------------|
| 20:1 | RPM    | 15 ~ 2,000 rpm |
|      | Torque | 5 ~ 70 N·cm    |



# IMPLANT LINE-UP



## Standard Implant Engine



SET UP

**COMPACT** **ECONOMIC** **DURABLE**

- 9 programmable memories
- Automatic overload protection
- Real time torque & RPM display
- Motor auto-calibration
- Multifunctional pedal control

|      |        |                |
|------|--------|----------------|
| 20:1 | RPM    | 30 ~ 2,500 rpm |
|      | Torque | 5 ~ 55 N·cm    |



## Hands-On Implant Engine



SET UP

**ECONOMIC** **BASIC** **HANDS-ON**

- Optimized for hands-on & seminar
- Sufficient torque up to 55 N·cm
- 3 programmable memories
- Powerful external pump system

|      |        |               |
|------|--------|---------------|
| 20:1 | RPM    | 30~ 1,500 rpm |
|      | Torque | 10 ~ 55 N·cm  |

