


DENTISSQ | SHORT IMPLANT



Safe and strong fixation

SQ Short Implant



Machined Surface

Strengthening the Bonding force
with soft tissue

Single Thread

Stable depth control

4 Blade Chip Pocket

Maximize Self-Tapping
Capabilities


20°

SQ Short Fixture Line-Up



Fixture Diameter	Length	Machine Collar	Code No.
5.0	6.0mm	2.0mm	DSSFR5004S
		1.0mm	DSSFR5005S
		0mm	DSSFR5006S


Fixture Package



Drilling Sequence


Place a fixture 0.5mm to 1mm deep in the lower part of the bone level.

Ø5.0 x 6mm




Normal bone

Hard bone



SQ SHORT&WIDE KIT



Tel. +82-1899-2804 | Fax. +82-53-583-2806
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea

©2020 DENTIS Co., Ltd.
All Right Reserved. Specifications are subject to change without notice.
Trademarks are the property of DENTIS Co., Ltd. Or their respective owners.
D-L-SQ-V1-202003-ENG



Tel. +1-323-677-4363 | Fax. +1-323-677-4366
6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +021-5111-3828 | Fax. +021-5111-3828
上海市长宁区中山西路933号2205室